

2016
Fourth Quarter
Issue 11
Volume 39

Shofu Solutions

Smart Products for a Healthy Smile



SHOFU

Visit www.shofu.com
or call 800.827.4638

is Issue

NEW! The Artistic Composite Kit

Developed by Frank J. Milnar,
DDS, AAACD

NEW! The Attachment Removal Kit for Clear Aligners

Consistent and Efficient
Dental Photography—
Stress-free!
NEW EyeSpecial C-II ebook

Events
January, February & March

Your Shofu Contacts

First Quarter Specials

Shofu Dental Corporation • San Marcos, CA

NEW! The Artistic Composite Kit

Designed for artisan's finishing of anterior composite restorations, the Artistic Composite Kit features a selection of Shofu's finest diamonds, carbides and stones. The instruments included in this kit are concomitant with "The Essence of Composite Art," an artistic restoration concept developed by Frank J. Milnar, DDS, AAACD. Arranged in a sequence of logical steps, these carefully-selected abrasives will help clinicians create natural-like aesthetics in an efficient manner.

The Artistic Composite Kit finishing technique supports minimally invasive dentistry. The protocol begins with creating a unified surface using fine and super fine Robot Diamonds in various shapes (PN 840F-1, 840V-1, 835F-1) which helps effectively prep the surface for further characterization.

In the following steps, with a fine-grit T&F Hybrid Point diamond (PN 0933-1) and egg-shaped Robot Carbide Finisher (PN 0718), three facial planes – the emergence profile, face of the tooth, and the incisal third — and secondary anatomy can be achieved. In areas of the restoration where light-dissipating architecture is needed, a flame-shaped Dura-Green Stone (PN 0060) will deliver predictable results. Other components employed for the final finishing step – the creation of macro and micro aesthetics, are the 8-, 16-, and 30- suted

Robot Carbide Finishers (PN 0770-3, 0771-3, 0772-3).

A finishing system with artisan's selected abrasives is an indispensable element in achieving natural aesthetics. Shofu's Artistic Composite Kit includes finest diamonds, carbides and stones, all designed to help create naturally appearing smile in a swift, predictable manner.

The Artistic Composite Kit

- Comprehensive finishing system inclusive of finest diamonds, carbides and stones
- Excellent for anterior composite restorations
- Supports minimally invasive dentistry

Suggested retail price: \$165

