

## Genie® Ultra Hydrophilic Impression Material

Performance without the high price tag!

- Ability to flow under pressure which helps capture detail beyond the prepared margin
- Your choice! Two set times – rapid 2:10 and standard 4:30 – in five viscosities

**BUY 6 boxes (2-50mL cartridges/box),  
GET 1 free!** (SNP16008)



## Genie® Bite Bite Registration Material

- Sets to a shore hardness of 90+
- Fluffy, moussy consistency stays put, without slumping
- Your choice! Superfast :60 and Fast :90



**BUY 3 boxes, GET 1 free!** (SNP16007)

## DuraShield® CV

5% Sodium Fluoride Clear Varnish



- Invisible<sup>1</sup>
- Significantly less wear time: 2 hours vs. 4-6 hours for other varnishes.<sup>1</sup>
- Thin, non-drip formula applies smoothly<sup>1</sup>
- Unique unit dose<sup>2</sup> designed for comfort and flexibility during application

**BUY 2 boxes (200-ct), GET 2 boxes (50-ct) free!**

(SNP16001)

**BUY 3 boxes (50-ct), GET 1 free!** (SNP16002)

## Upgrade®

Disposable Prophy Angles

- Great flare each time for better coverage and to aid in cleaning
- Webbed and ribbed prophy cup adapts to tooth contours and provides very good splatter control
- Lubricated gearing for smooth performance
- Multi-color pack

**BUY 4 packages (100-ct), GET 1 free!** (SNP16005)

**BUY 1 package (1200-ct), GET 3 packages (100-ct) free!** (SNP16004)



## 3-Way®

Disposable Impression Trays



**BUY 3 boxes, GET 1 free!**

(SNP16006)

## Silgimix™

Alginate Replacement Impression Material

Ideal for making impression study models, orthodontic models, counter impressions, provisionals, bleaching trays and more.

- Pour Now, Pour Later, Pour Again! Impressions can be poured multiple times without loss of accuracy.
- Ideal for making impressions that will be scanned as part of a CAD/CAM process.



## Versa-Temp® 2

Temporary Crown and Bridge Resin



Combines superb performance and natural-looking esthetics.

- Up to 70% stronger<sup>3</sup> than leading competitors – 450 MPa compressive strength
- 4 minute total set time
- Natural Fluorescence
- Easy to polish with an alcohol wipe (e.g. Isopropyl)

**BUY 3 single packs, GET 1 free!** (SNP16009)

Excludes Bulk Kits

## Topex® Prophy Paste

An assortment of great tasting flavors sure to satisfy any patient.



**BUY 4 boxes (200-ct), GET 1 free!** (SNP16003)

## Moist SURE™ Soothing Hand Hygiene

A complete line of hand hygiene products. Sanitizers, soap, lotion and a dispenser—all designed to care for professional hands and encourage compliance.

**BUY 1 Auto Dispenser and 1 refill (33.8oz), GET 1 refill (33.8oz) free!** (SNP16014)

**BUY 4 pump bottles (15oz) or sanitizers (9oz), GET 1 of the same free!** (SNP16015)



## Pro-Sonic® Ultrasonic Cleaners

Efficient, effective and economical—and come in a variety of sizes to suit any office. Backed by one of the best warranties in the industry!



**BUY 1, GET 1 ReSURGE bottle (33.8oz) free!**

(Shipped with order)

## ReSURGE® Instrument Cleaning Solution

- Unique formula helps protect your instruments' passive layer (outer coating), reducing the potential for stains, spots and even rust.<sup>1</sup>
- Biodegradable



**BUY 1 bottle (33.8 oz), GET 1 ASSURE Plus™ Sterilization Pouches 3.5 x 10 (200-ct) free!** (SNP16019)

## ASSURE Plus™

Self-Seal Sterilization Pouches



The only sterilization pouch with a built in peel-off internal indicator, as recommended by the CDC.

- Unique peel-off indicator
- Dual external heat and gas sterilization indicators
- Strong, medical grade paper that allows steam and gas penetration
- Easy to open overlap design

**BUY 5 boxes, GET 1 free!** (SNP16010)

## Sporox II®

Sterilizing and Disinfecting Solution



Powerful hydrogen peroxide-based sterilant and high-level disinfectant for heat-sensitive instruments.

- Ready to use - no mixing, heating or activation required
- No noxious odors, safe for your staff
- Will not bond infectious proteins to instruments like glutaraldehydes
- Sterilization is achieved in 6 hours, high-level disinfection in 30 minutes
- Solution can be re-used for up to 21 days

**BUY 3 gallons, GET 1 ProSoak Tray free!** (SNP16020)

## FlashTips®

Disposable Air/Water Syringe Tips



Designed to prevent cross contamination

- No converters needed on most syringes
- Locking indicator that assures tip is in place
- Built strong for cheek retraction

**BUY 1 bag (1200-ct), GET 1 (250-ct) free!** (SNP16013)

**BUY 4 bags (250-ct), GET 1 free!** (SNP16012)

## Com-Fit Plush™ Masks

Com-Fit Plush™ masks contain an innovative inner material that feels silky soft. They're hypoallergenic for sensitive skin and available in three levels of fluid resistance.<sup>2</sup>



**BUY any 4 boxes of Com-Fit Masks, GET 1 box of Com-Fit Plush free!** (SNP16011)

## Purevac®

Evacuation System Cleaners



Outperforms leading competitors in percentage of soil removal in a cleaning study when using a Purevac® solution and the Purevac® Autowalkabout Dispenser. See the study at purevac.net.

**BUY 2 bottles (2-liter), GET 2 boxes of Com-Fit Plush Masks free!** (SNP15072)

**BUY 2 bottles (5-liter), GET 4 boxes of Com-Fit Plush Masks free!** (SNP15073)

## VoloWipes™ Disinfectant Wipes

Sized to fit your hand, and the job, perfectly!

- 18.5% larger than standard wipes
- Fast disinfection against microorganisms! Tuberculocidal, virucidal, and bactericidal in just two minutes



**BUY 12, GET 2 free!** (SNP16021)

1. Instrument protection study 2012 2. Tested per ASTM F1862.

Purchases must be made between 1/1/16 and 3/31/16. Purchases must be on one invoice. Offers must be redeemed by 4/30/16. To redeem your FREE GOODS, mail your invoice to: Sultan Healthcare, Attn: FREE GOODS, 1301 Smile Way, York, PA 17404-0807, or fax to: Sultan Healthcare, Attn: FREE GOODS, 1-201-968-1409. Allow 4-6 weeks for delivery. Offers valid in the 50 United States and DC only. Limit 5 redemptions per doctor. Offers may not be combined with any other Sultan offers or contract agreements. Data on file. ©2015 Sultan Healthcare. All rights reserved.

