

# Glaze Ceramics Like Never Before!

Buy any 2 like items of Enamelite® Products, Get 1 of equal or lesser value at no charge!



## Enamelite® Ceramic Spray

- Accurate, evenly-sprayed coat
- Glaze single or multiple cases in seconds
- Waste is reduced to a minimum

CAD/CAM Scanning Spray (8g/3.1 oz) #5100176

Low Fusing Small (56g/2.0oz) #5100180


Low Fusing Medium (111g/3.9oz) #5100179

Low Fusing Fluorescent Small (56g/2oz.) #5100181

# The Leading Laminates in the Dental Industry

\*Please visit [dental.keystoneindustries.com](http://dental.keystoneindustries.com) for item numbers

Buy any 2 boxes of Pro-form® or EVA, Get 1 box of equal or lesser value at no charge!



## EVA Resin Sheets

- Comes in thicknesses of: .160 and .120 (youth mouthguard)
- Available in 5"x 5" (127 x 127mm) squares
- Comes in colors of: clear, red, blue, green, yellow, black, orange, white
- Available in packs of 25

## Soft Clear EVA

- Used for the fabrication of customer fluoride trays, bleaching trays and model duplication
- Available in 5"x 5" (127 x 127mm) squares
- Comes in thicknesses of: 0.040, 0.060, 0.080, 0.120, 0.160
- Available in packs of 25

## Pro-form® Laminates

- Available in 3 sizes:
  - 5"x 5" (127 x 127mm) squares
  - 125mm, 120mm round
- Pro-form Laminate & Fun Guards are available in packs of 12
- Tie-Dye & Glitter Guards are available in packs of 6
- Assorted color packs available



This offer cannot be combined with any other promotion or discount. Other restrictions may apply. Customer will be invoiced for purchased goods from preferred distributor and no charge goods will ship with the same order. Offer valid from 10/01/2018—01/04/2010. Valid in the US only. Not responsible for typographical errors.